Calibrador β_2 -Microglobulina

MonlabTest®

Turbidimetría

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Turbidimetry

Quantitative determination of β_2 -Microglobulin (β_2 -m)

Only for professional *in vitro* diagnostic use.
Store at 2-8°C.

PRODUCT CHARACTERISTICS

The calibrator is a lyophilised serum from human origin used to determinate β_2 -m concentration in human samples by turbidimetry or nephelometry methods.

REAGENTS

β_2-m Calibrator	Calibrator. Human serum. β_2 -m concentration is stated on the vial label.
--	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator is calibrated against the 1st International Standard ERM-DA470k/IFCC.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Wait for 10 minutes before use.

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 1 month at 2-8°C or 3 months at -20°C.

PROCEDURE

To be used in turbidimetric and nephelometric assays.

PRESENTATION

MO-165053

1 x 1 mL