Giardia+ MonlabTest®

MO-804008 20 TESTS

One Step Giardia Antigen Test Device

IVD

A rapid, one step test for the qualitative detection of *Giardia* (CWP1 and/or α 1-giardin) in human faeces.

For professional in vitro diagnostic use only.

INTENDED USE

Giardia+ MonlabTest® test is a rapid chromatographic immunoassay for the qualitative detection of *Giardia* (CWP1 and/or α 1-giardin) in human faeces specimens to aid in the diagnosis of giardiasis.

SYNTHESIS

Giardiasis is a diarrhoeal illness seen throughout the world. It is caused by a flagellate protozoan parasite, *Giardia intestinalis*, also known as *G. lamblia* and *G. duodenalis*.

Giardia is a common cause of gastrointestinal disturbance in both high and low-income countries. The incidence of *Giardia* is generally higher in low-income countries (e.g. many countries of Africa, Asia, and South and Central America) where access to clean water and basic sanitation is lacking. Nearly all children in this setting will acquire *Giardia* at some point in their childhood, and the prevalence of the parasite in young children can be as high as 10%-30%. In areas such as Western Europe and the United States of America, *Giardia* infection is associated with ingestion of contaminated water, person-to-person spread, recent foreign travel, and recreational swimming. *Giardia* may be cause of 2%-5% of cases of diarrhoea in high-income countries.

Giardia lamblia is a flagellated unicellular eukaryotic microorganism that commonly causes diarrheal disease worldwide. It is the most common cause of waterborne diarrheal outbreaks in the United States and is occasionally considered a cause of foodborne diarrhea.

In developing countries, there is a very high prevalence and incidence of infection, and data suggest that chronic giardiasis can cause long-term stunting. *Giardia* species have important stages in the life cycle. A host infection begins when the cyst is swallowed with contaminated water. After exposure to the acidic environment of the stomach, the cyst is excreted in trophozoites in the proximal small intestine. Trophozoite is the vegetative form and replicates in the small intestine, where it causes symptoms of diarrhea and malabsorption.

PRINCIPLE

The Giardia MonlabTest® is a qualitative immunoassay for the detection of *Giardia* (CWP1 and/or α 1-giardin) in human faeces samples. The membrane is pre-coated with antibodies against *Giardia* (CWP1 and/or α 1-giardin) on the test line region. During testing, the sample reacts with the particle coated with anti-*Giardia* (CWP1 and/or α 1-giardin) antibodies which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. In the case of a positive result the specific antibodies present on the membrane will react with the mixture conjugates and generate coloured lines. A green coloured band always appears in the control line and serves as verification that sufficient volume was added that proper flow was obtained and as an internal control for the reagents.

PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pack until use.
- Do not use the test if pack is damaged.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test should be discarded in a proper biohazard container after testing.

- The test must be carried out within 2 hours of opening the sealed bag.
- The presence of yellowed lines in the results window (zone of control and test lines) that are visible before using the test are completely normal. These lines do not mean test functionality failure.

STORAGE AND STABILITY

Store as packaged in the sealed pack either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pack. The test must remain in the sealed pack until use. Do not freeze.

MATERIAL PROVIDED

- 20 Tests
- Instructions for use
- 20 specimen collection vial with buffer

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

SPECIMEN COLLECTION AND PREPARATION

Collect sufficient quantity of faeces (1-2 g or mL for liquid sample). Stool samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator 2-4°C/36-40°F) for 1-2 days prior to testing. For longer storage the specimen must be kept frozen at -20°C/-4°F. In this case, the sample will be totally thawed, and brought to room temperature before testing.

PROCEDURED

To process the collected stool samples (see illustration 1):

Use a separate specimen collection vial for each sample with 1 mL of the buffer. Unscrew the cap of the vial and introduce the stick four times and cover the propeller of the stick to take enough sample, trying not to exceed the area of the helix and thus avoid erroneous results.

Close the vial with the buffer and stool sample. Shake the vial in order to assure good sample dispersion. For liquid stool samples, aspirate the faecal specimen with a dropper and add 125µL into the specimen collection vial with buffer.

Test Procedure (see illustration 2)

Allow the tests, stool samples and buffer to reach to room temperature (15-30°C/59-86°F) prior to testing. Do not open pouches until ready to perform the assay.

1. Remove the Giardia+ MonlabTest® from its sealed pouch and use it as soon as possible.
2. Shake the specimen collection vial to assure good sample dispersion. Break off the cap of the vial.
3. Use a separate device for each sample. Dispense exactly 3 drops into the specimen well (S). Start the timer.
4. Read the result at 10 minutes after dispensing the sample.

Illustration 1



Illustration 2

INTERPRETATION OF RESULTS

Illustration 3

POSITIVE: Two lines appears across the central window. In the result line region, a red test line called Giardia test line, marked in the illustration 3 with the letter T, and in the control line region, a green control line marked in the illustration 3 with the letter C.

NEGATIVE: Only one green band appears across the control line region marked with the letter C at the illustration 3 (control line).

INVALID: A total absence of the green control coloured band regardless the appearance or not of the red test line.

Note: Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (I) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test:
- A green line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- Giardia+ MonlabTest® will only indicate the presence of parasites in the specimen (qualitative detection) and should be used for the detection of *Giardia* (CWPI and/or $\alpha 1$ -giardin) in faeces specimens only. Neither the quantitative value nor the rate of increase in antigen concentration can be determined by this test.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- Do not use specimens treated with solutions containing formaldehyde or its derivatives.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of giardiasis.
- After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- This test provides a presumptive diagnosis of giardiasis. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- Mucous and/or bloody stool samples may cause non-specific reactions in the test. Mucous and/or bloody stool samples that are positive should be confirmed by other techniques.

EXPECTED VALUES

Giardia is prevalent throughout the world, including temperate, high-income countries, such as the UK and the United States. Several studies have examined acquisition of giardiasis in international travellers.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

It was performed an evaluation using Giardia+ MonlabTest+ vs Evaluation Criteria.

The Evaluation Criteria is: Giardia+ MonlabTest® is compared with other two LF tests in the market and when the results obtained with any of the test are not concordant, the samples must be evaluated with a qPCR test of the market (*Cryptosporidium*, giardia & *E. histolytica* Real Time Detection kit).

IC test: Giardia+ MonlabTest®	Evaluation criteria			Total
	+	-		
+	44	2	46	
-	1	79	80	
Total	45	81	126	

Giardia+ MonlabTest® vs Evaluation criteria		
Sens.	97.8%	88.2 - 99.9%
Esp.	97.5%	91.4 - 99.7%
VPP	95.7%	85.2 - 99.5%
VPN	98.8%	93.2 - 100.0%

Cross-Reactivity and Interferences

It was performed an evaluation to determine the cross reactivity and interferences of Giardia+ MonlabTest®. There are no cross-reactions with gastrointestinal pathogens, other organisms, substances and/or fecal markers occasionally present in the stool.

<i>Adenovirus</i>	<i>Lactoferrina humana</i>
<i>Astrovirus</i>	<i>Legionella</i>
<i>Calprotectin</i>	<i>Listeria monocytogenes</i>
<i>Campylobacter jejuni</i>	<i>Norovirus GI/GII</i>
<i>Clostridium difficile</i> GDH/ <i>Toxina A/Toxina B</i>	<i>Rotavirus</i>
<i>Clostridium perfringens</i>	<i>Salmonella enteritidis/ paratyphi</i> <i>/typhi/typhimurium</i>
<i>Cryptosporidium paroum</i>	<i>Shigella boydii</i> <i>/dysenteriae/flexneri/sonnei</i>
<i>E. coli</i> O157 / O26 / O111	<i>Streptococcus</i> <i>pyogenes/Pneumococcal</i>
<i>Entamoeba histolytica</i>	<i>Transferrina humana</i>
<i>H. pylori</i>	<i>Yersinia enterocolitica</i> O:3/O:9
<i>Hemoglobina</i> <i>humana/bovina/cerdo</i>	

REFERENCES

- Hill DR, Nash TE. Intestinal Flagellate and Ciliate Infections. In: Guerrant RL, Walker DH, Weller PF, eds. Tropical Infectious Diseases. Principles, Pathogens & Practice. 2nd ed. Elsevier, Philadelphia. 2006:984-8.
- Copue S, Delabre K, Pouillot R et al. Detection of *Cryptosporidium*, *Giardia* and *Enterocytozoombieneusi* in surface water, including recreational areas: a one year prospective study. FEMS Immunol Med Microbiol. 2006; 47:351-9.
- Stuart JM, Orr HJ, Warburton FG, et al. Risk Factors for Sporadic Giardiasis: A Case-Control Study in Southwestern England. Emerg. Infect Dis. 2003; 9, 2.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> test		Keep dry
	Catalogue Code		Temperature limitation
	Lot number		Use by

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1. Hill DR, Nash TE. Intestinal Flagellate and Ciliate Infections. In: Guerrant RL, Walker DH, Weller PF, eds. Tropical Infectious Diseases. Principles, Pathogens & Practice. 2nd ed. Elsevier, Philadelphia. 2006:984-8.
2. Copue S, Delabre K, Pouillot R et al. Detection of *Cryptosporidium*, *Giardia* and *Enterocytozoon bieneusi* in surface water, including recreational areas: a one year prospective study: FEMS Immunol Med Microbiol. 2006; 47:351-9.
3. Stuart JM, Orr HJ, Warburton FG, et al. Risk Factors for Sporadic Giardiasis: A Case-Control Study in Southwestern England. Emerg Infect Dis. 2003; 9, 2.

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